

Anti-REST antibody (760-1060) (STJ116170)

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-REST (760-1060) is suitable for use in Immunofluorescence.

Applications IF Host/Source Rabbit Reactivity Human

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID

Concentration

Conjugation Unconjugated Purification Affinity purification Dilution Range IF 1:50-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 5978 Gene Symbol REST

Uniprot ID REST_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 760-1060 of human REST/NRSF

(NP_001180437.1).

Immunogen Region 760-1060

Specificity Immunogen Sequence

