

Anti-EIF5B antibody (901-1220) (STJ116285)

STJ116285

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EIF5B (901-1220) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:10-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9669
Gene Symbol	EIF5B
Uniprot ID	IF2P_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 901-1220 of human EIF5B (NP_056988.3).
Immunogen Region	901-1220
Specificity	
Immunogen Sequence	