

Anti-TRIM13 antibody (40-310) (STJ116298)

STJ116298

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TRIM13 (40-310) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:20-1:50
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10206
Gene Symbol	TRIM13
Uniprot ID	TRI13_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-310 of human TRIM13 (NP_001007279.1).
Immunogen Region	40-310
Specificity	
Immunogen Sequence	