

Anti-TPD52L1 antibody (1-144) (STJ116303)

STJ116303

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TPD52L1 (1-144) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:20-1:50
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7164
Gene Symbol	TPD52L1
Uniprot ID	TPD53_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human TPD52L1 (NP_001003396.1).
Immunogen Region	1-144
Specificity	
Immunogen Sequence	