

Anti-SPAG5 antibody (1-300) (STJ116304)

STJ116304

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SPAG5 (1-300) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:20-1:50
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10615
Gene Symbol	SPAG5
Uniprot ID	SPAG5_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SPAG5 (NP_006452.3).
Immunogen Region	1-300
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of HeLa cells using SPAG5 antibody (STJ116304). Blue: DAPI for nuclear staining.