

## Anti-DDIT3 antibody (STJ116305)

## **GENERAL INFORMATION**

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-DDIT3/CHOP is suitable for use in Western Blot, Immunohistochemistry and

Immunofluorescence.

Applications WB, IHC, IF Host/Source Rabbit

Reactivity Human, Mouse, Rat

## **PRODUCT PROPERTIES**

Clonality Polyclonal

Clone ID Concentration

Conjugation Unconjugated Purification Affinity purification Dilution Range WB 1:500-1:2000

IF 1:50-1:200 IHC 1:50-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

## TARGET INFORMATION

Gene ID 1649 Gene Symbol DDIT3

Uniprot ID DDIT3\_HUMAN

Immunogen Recombinant protein of human DDIT3/CHOP Immunogen Region Specificity Immunogen Sequence







