

Anti-CDCA4 antibody (1-241) (STJ116328)

STJ116328

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CDCA4 (1-241) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55038
Gene Symbol	CDCA4
Uniprot ID	CDCA4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-241 of human CDCA4 (NP_060425.2).
Immunogen Region	1-241
Specificity	
Immunogen Sequence	