

Anti-MCTS1 antibody (1-181) (STJ116340)

STJ116340

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MCTS1 (1-181) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	28985
Gene Symbol	MCTS1
Uniprot ID	MCTS1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human MCTS1 (NP_054779.1).
Immunogen Region	1-181
Specificity	
Immunogen Sequence	