

Anti-AMY2B antibody (16-320) (STJ116438)

STJ116438

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AMY2B (16-320) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	280
Gene Symbol	AMY2B
Uniprot ID	AMY2B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-320 of human AMY2B (NP_066188.1).
Immunogen Region	16-320
Specificity	
Immunogen Sequence	