

Anti-AMZ2 antibody (81-360) (STJ116545)

STJ116545

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AMZ2 (81-360) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51321
Gene Symbol	AMZ2
Uniprot ID	AMZ2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 81-360 of human AMZ2 (NP_057711.3).
Immunogen Region	81-360
Specificity	
Immunogen Sequence	