

Anti-ST7 antibody (200-290) (STJ116596)

STJ116596

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ST7 (200-290) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7982
Gene Symbol	ST7
Uniprot ID	ST7_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-290 of human ST7 (NP_068708.1).
Immunogen Region	200-290
Specificity	
Immunogen Sequence	