

Anti-TRIM2 antibody (100-300) (STJ116606)

STJ116606

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TRIM2 (100-300) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23321
Gene Symbol	TRIM2
Uniprot ID	TRIM2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human TRIM2 (NP_056086.2).
Immunogen Region	100-300
Specificity	
Immunogen Sequence	