

Anti-RPUSD2 antibody (286-545) (STJ116609)

STJ116609

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RPUSD2 (286-545) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27079
Gene Symbol	RPUSD2
Uniprot ID	RUSD2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 286-545 of human RPUSD2 (NP_689473.1).
Immunogen Region	286-545
Specificity	
Immunogen Sequence	