

Anti-ZFP82 antibody (50-170) (STJ116671)

STJ116671

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZFP82 (50-170) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	284406
Gene Symbol	ZFP82
Uniprot ID	ZFP82_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-170 of human ZFP82 (NP_597723.1).
Immunogen Region	50-170
Specificity	
Immunogen Sequence	