

Anti-CLDN8 antibody (29-81) (STJ116680)

STJ116680

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CLDN8 (29-81) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9073
Gene Symbol	CLDN8
Uniprot ID	CLD8_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-81 of human CLDN8 (NP_955360.1).
Immunogen Region	29-81
Specificity	
Immunogen Sequence	