

Anti-ASF1A antibody (100-204) (STJ116719)

STJ116719

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ASF1A (100-204) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25842
Gene Symbol	ASF1A
Uniprot ID	ASF1A_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-204 of human ASF1A (NP_054753.1).
Immunogen Region	100-204
Specificity	
Immunogen Sequence	