

Anti-CDSN antibody (100-200) (STJ116779)

STJ116779

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CDSN (100-200) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1041
Gene Symbol	CDSN
Uniprot ID	CDSN_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CDSN (NP_001255.3).
Immunogen Region	100-200
Specificity	
Immunogen Sequence	