

Anti-CDH18 antibody (300-500) (STJ116851)

STJ116851

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CDH18 (300-500) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1016
Gene Symbol	CDH18
Uniprot ID	CAD18_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-500 of human CDH18 (NP_004925.1).
Immunogen Region	300-500
Specificity	
Immunogen Sequence	