

Anti-PDGFD antibody (185-364) (STJ116859)

STJ116859

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDGFD (185-364) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80310
Gene Symbol	PDGFD
Uniprot ID	PDGFD_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 185-364 of human PDGFD (NP_149126.1).
Immunogen Region	185-364
Specificity	
Immunogen Sequence	