

Anti-PFKFB3 antibody (1-210) (STJ116964)

STJ116964

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PFKFB3 (1-210) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5209
Gene Symbol	PFKFB3
Uniprot ID	F263_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human PFKFB3 (NP_004557.1).
Immunogen Region	1-210
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of U2OS cells using PFKFB3 antibody (STJ116964). Blue: DAPI for nuclear staining.