

Anti-ZNF141 antibody (405-474) (STJ116999)

STJ116999

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZNF141 (405-474) is suitable for use in Immunohistochemistry and Immunofluorescence.
Applications	IHC, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7700
Gene Symbol	ZNF141
Uniprot ID	ZN141_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 405-474 of human ZNF141 (NP_003432.1).
Immunogen Region	405-474
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded mouse lung using ZNF141 antibody (STJ116999) at dilution of 1:100 (40x lens).