

Anti-PLA2G6 antibody (350-450) (STJ117008)

STJ117008

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PLA2G6 (350-450) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53357
Gene Symbol	Pla2g6
Uniprot ID	PLPL9_MOUSE
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of mouse PLA2G6 (NP_001185952.1).
Immunogen Region	350-450
Specificity	
Immunogen Sequence	