

Anti-LRP10 antibody (17-120) (STJ117067)

STJ117067

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LRP10 (17-120) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26020
Gene Symbol	LRP10
Uniprot ID	LRP10_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-120 of human LRP10 (NP_054764.2).
Immunogen Region	17-120
Specificity	
Immunogen Sequence	