

Anti-SLC40A1 antibody (200-300) (STJ117085)

STJ117085

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC40A1 (200-300) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	30061
Gene Symbol	SLC40A1
Uniprot ID	S40A1_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human SLC40A1 (NP_055400.1).
Immunogen Region	200-300
Specificity	
Immunogen Sequence	