

Anti-DDX59 antibody (1-250) (STJ117139)

STJ117139

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DDX59 (1-250) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83479
Gene Symbol	DDX59
Uniprot ID	DDX59_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DDX59 (NP_001026895.2).
Immunogen Region	1-250
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded rat testis using DDX59 antibody (STJ117139) at dilution of 1:100 (40x lens).