

Anti-AMICA1 antibody (285-384) (STJ117157)

STJ117157

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-AMICA1 (285-384) is suitable for use in Western Blot. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB 1:500-1:2000 |
| Formulation | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. |
| Isotype | IgG |
| Storage Instruction | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 120425 |
| Gene Symbol | JAML |
| Uniprot ID | JAML_HUMAN |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 285-384 of human AMICA1 (NP_694938.2). |
| Immunogen Region | 285-384 |
| Specificity | |
| Immunogen Sequence | |