

Anti-ACP1 antibody (15-144) (STJ117214)

STJ117214

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ACP1 (15-144) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	52
Gene Symbol	ACP1
Uniprot ID	PPAC_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-144 of human ACP1 (NP_004291.1).
Immunogen Region	15-144
Specificity	
Immunogen Sequence	