

Anti-ADAM8 antibody (17-300) (STJ117216)

STJ117216

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ADAM8 (17-300) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	101
Gene Symbol	ADAM8
Uniprot ID	ADAM8_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-300 of human ADAM8 (NP_001100.3).
Immunogen Region	17-300
Specificity	
Immunogen Sequence	