

Anti-BNIP3 antibody (1-140) (STJ117226)

STJ117226

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BNIP3 (1-140) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	664
Gene Symbol	BNIP3
Uniprot ID	BNIP3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human BNIP3 (NP_004043.3).
Immunogen Region	1-140
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of U2OS cells using BNIP3 antibody (STJ117226) at dilution of 1:100. Blue: DAPI for nuclear staining.