

Anti-FUT3 antibody (1-100) (STJ117252)

STJ117252

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-FUT3 (1-100) is suitable for use in Western Blot and Immunofluorescence.

Applications WB, IF Host/Source Rabbit Reactivity Human, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID Concentration

Conjugation
Purification Unconjugated
Affinity purification
Dilution Range WB 1:500-1:2000

IF 1:50-1:100

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG


Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 2525 Gene Symbol FUT3 Uniprot ID FUT3_HUMAN

Immunogen A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FUT3 (NP_000140.1).

Immunogen Region 1-100 Specificity Immunogen Sequence



Western blot analysis of extracts of various cell lines, using FUT3 antibody (STJ117252) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 39% nonfat dry milk in TBST. Detection: