

Anti-NPBWR1 antibody (150-250) (STJ117255)

STJ117255

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NPBWR1 (150-250) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2831
Gene Symbol	NPBWR1
Uniprot ID	NPBW1_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human NPBWR1 (NP_005276.2).
Immunogen Region	150-250
Specificity	
Immunogen Sequence	