

Anti-IFNW1 antibody (22-195) (STJ117261)

STJ117261

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IFNW1 (22-195) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3467
Gene Symbol	IFNW1
Uniprot ID	IFNW1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-195 of human IFNW1 (NP_002168.1).
Immunogen Region	22-195
Specificity	
Immunogen Sequence	