

Anti-NPY antibody (1-97) (STJ117270)

STJ117270

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NPY (1-97) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4852
Gene Symbol	NPY
Uniprot ID	NPY_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-97 of human NPY (NP_000896.1).
Immunogen Region	1-97
Specificity	
Immunogen Sequence	