

Anti-PSG2 antibody (200-300) (STJ117280)

STJ117280

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PSG2 (200-300) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5670
Gene Symbol	PSG2
Uniprot ID	PSG2_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PSG2 (NP_112536.2).
Immunogen Region	200-300
Specificity	
Immunogen Sequence	