

Anti-CXCR4 antibody (1-100) (STJ117306)

STJ117306

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CXCR4 (1-100) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7852
Gene Symbol	CXCR4
Uniprot ID	CXCR4_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CXCR4 (NP_003458.1).
Immunogen Region	1-100
Specificity	
Immunogen Sequence	