

Anti-RABGAP1L antibody (50-210) (STJ117318)

STJ117318

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RABGAP1L (50-210) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9910
Gene Symbol	RABGAP1L
Uniprot ID	RBG10_HUMAN RBG1L_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-210 of human RABGAP1L (NP_055672.3).
Immunogen Region	50-210
Specificity	
Immunogen Sequence	