

Anti-CEP170 antibody (1505-1584) (STJ117565)

STJ117565

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CEP170 (1505-1584) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9859
Gene Symbol	CEP170
Uniprot ID	CE170_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1505-1584 of human CEP170 (NP_055627.2).
Immunogen Region	1505-1584
Specificity	
Immunogen Sequence	