

Anti-ARFGF2 antibody (220-370) (STJ117583)

STJ117583

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ARFGF2 (220-370) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10564
Gene Symbol	ARFGF2
Uniprot ID	BIG2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-370 of human ARFGF2 (NP_006411.2).
Immunogen Region	220-370
Specificity	
Immunogen Sequence	