

Anti-FNDC3A antibody (900-1170) (STJ117599)

STJ117599

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FNDC3A (900-1170) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22862
Gene Symbol	FNDC3A
Uniprot ID	FND3A_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1170 of human FNDC3A (NP_001265367.1).
Immunogen Region	900-1170
Specificity	
Immunogen Sequence	