

Anti-KIAA1324 antibody (42-300) (STJ117685)

STJ117685

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KIAA1324 (42-300) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57535
Gene Symbol	ELAPOR1
Uniprot ID	ELAP1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-300 of human KIAA1324 (NP_001253977.1).
Immunogen Region	42-300
Specificity	
Immunogen Sequence	