

Anti-MED28 antibody (1-178) (STJ117712)

STJ117712

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MED28 (1-178) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80306
Gene Symbol	MED28
Uniprot ID	MED28_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human MED28 (NP_079481.2).
Immunogen Region	1-178
Specificity	
Immunogen Sequence	