

Anti-FNIP1 antibody (270-350) (STJ117747)

STJ117747

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FNIP1 (270-350) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	96459
Gene Symbol	FNIP1
Uniprot ID	FNIP1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-350 of human FNIP1 (NP_588613.2).
Immunogen Region	270-350
Specificity	
Immunogen Sequence	