

Anti-IFNA4 antibody (24-189) (STJ117849)

STJ117849

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IFNA4 (24-189) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3441
Gene Symbol	IFNA4
Uniprot ID	IFNA4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-189 of human IFNA4 (NP_066546.1).
Immunogen Region	24-189
Specificity	
Immunogen Sequence	