

Anti-C1QC antibody (50-150) (STJ117864)

STJ117864

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C1QC (50-150) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	714
Gene Symbol	C1QC
Uniprot ID	C1QC_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human C1QC (NP_001107573.1).
Immunogen Region	50-150
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded human colon using C1QC antibody (STJ117864) at dilution of 1:100 (40x lens).