

Anti-Phospho-PRKCA-S319 antibody (STJ117877)

STJ117877

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-Phospho-PKCalpha-S319 is suitable for use in Western Blot.

Applications WB Host/Source Rabbit

Reactivity Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID

Concentration

Conjugation Unconjugated Purification Affinity purification Dilution Range WB 1:500-1:2000

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG


Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 5578
Gene Symbol PRKCA
Uniprot ID KPCA_HUMAN

Immunogen
Immunogen Region
Specificity
Immunogen Sequence

Immunogen A synthetic phosphorylated peptide around S319 of human PKC alpha (NP_002728.2).



Western blot analysis of extracts of Hela cells, using Phospho-PKC alpha-S319 antibody (STJ17877) at 1:1000 dilution. Hela cells were treated by PMM/TPA (200 nM) at 37° Cfor 30 minutes after serum-starvatior overnight. Hela cells were treated by ATP (5 mM) at 30° Cfor 1 hour. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins 250up per lane. Blocking buffer: 3% honfatt dry milk ir