

## Anti-Phospho-LIMK1/LIMK2-T505 antibody (STJ117878)

STJ117878

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-LIMK1-T508/LIMK2-T505 is suitable for use in Western Blot.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB 1:500-1:2000
<b>Formulation</b>	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store in a freezer at -20°C and avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3984</a> <a href="#">3985</a>
<b>Gene Symbol</b>	<a href="#">LIMK1</a> <a href="#">LIMK2</a>
<b>Uniprot ID</b>	<a href="#">LIMK1_HUMAN</a> <a href="#">LIMK2_HUMAN</a>
<b>Immunogen</b>	A synthetic phosphorylated peptide around T508 of human LIMK1 (NP_002305.1).
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	