

Anti-COL11A1 antibody (260-430) (STJ118005)

STJ118005

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-COL11A1 (260-430) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1301
Gene Symbol	COL11A1
Uniprot ID	COBA1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-430 of human COL11A1 (NP_542196.2).
Immunogen Region	260-430
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded rat heart using COL11A1 antibody (STJ118005) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse heart using COL11A1 antibody (STJ118005) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081