

Anti-ACADM antibody (26-421) (STJ118030)

STJ118030

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ACADM (26-421) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	34
Gene Symbol	ACADM
Uniprot ID	ACADM_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-421 of human ACADM (NP_000007.1).
Immunogen Region	26-421
Specificity	
Immunogen Sequence	