

Anti-DICER1 antibody (350-450) (STJ118244)

STJ118244

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DICER1 (350-450) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23405
Gene Symbol	DICER1
Uniprot ID	DICER_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human DICER1 (NP_001258211.1).
Immunogen Region	350-450
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of U2OS cells using DICER1 antibody (STJ118244) at dilution of 1:100. Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081