

Anti-SIRT4 antibody (1-200) (STJ118245)

STJ118245

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SIRT4 (1-200) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23409
Gene Symbol	SIRT4
Uniprot ID	SIR4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human SIRT4 (NP_036372.1).
Immunogen Region	1-200
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of U2OS cells using SIRT4 antibody (STJ118245) at dilution of 1:100. Blue: DAPI for nuclear staining.